



# VIELIGHT MIP

MULTIMODAL INTRANASAL PHOTOBIOMODULATION



[WWW.VIELIGHT.COM](http://WWW.VIELIGHT.COM)

# Why the **Nose**?

## THE NOSE IS MORE THAN AN AIRWAY.

Its unique anatomy makes it a distinctive route for intranasal photobiomodulation.



### BRAIN UNDERSIDE

The roof of the nasal cavity sits directly beneath the brain, making it a uniquely positioned route for intranasal photobiomodulation.



### RICH VASCULAR NETWORK

The nasal cavity contains a dense network of blood vessels, providing a pathway for systemic photobiomodulation.



### AIR PASSAGEWAYS


As a primary entry point for inhaled air, the nose continually interacts with airborne particles, moisture, and microorganisms.



### EYES & OPTIC NERVES

The nasal cavity is closely connected to nearby ocular structures, including pathways associated with the eyes and optic nerves.

**A SMALL SPACE WITH  
POWERFUL CONNECTIONS**





# What Is the Vielight MIP?

The Vielight MIP combines carefully selected wavelengths of light to support brain and systemic wellness in one compact device.

With two intranasal applicators and a portable controller, MIP is designed to support:

- ✓ Gentle brain stimulation
- ✓ Immunity
- ✓ Oxygenation
- ✓ Systemic wellness

The 470 nm configuration also supports intranasal microbial sanitization.

## WHAT IS PHOTOBIO-MODULATION?

Photobiomodulation (PBM) uses specific wavelengths of light to support cellular activity.

Light is absorbed by mitochondria — the energy centers of cells — helping support ATP production and normal cellular performance.



# Four Wavelengths

## DIFFERENT WAVELENGTHS. DIFFERENT ROLES.

MIP uses four distinct wavelengths, each designed to support a different pathway and wellness function.




**470 nm**

**Blue LED**

### **Intranasal microbial sanitization**

Designed to support microbial sanitization within the nasal cavity.




**633 nm**

**Red LED**

### **Systemic stimulation**

Designed to support systemic wellness through interaction with the nasal vasculature.



**655 nm**

**Red Laser**

### **Focused systemic stimulation**

A focused red-light option designed to support systemic wellness through the nasal vasculature.



**810 nm**

**Near-Infrared LED**

### **Pulses at 10Hz or 40Hz**

Designed to support neural mitochondria and blood flow.

**10Hz (Alpha)** for relaxation and overall brain healing

**40Hz (Gamma)** for focus, memory, and cognition



**One platform.  
Four targeted roles.  
Four wavelengths.  
Four distinct pathways.  
One wellness system.**

# Find Your MIP

## CHOOSE THE CONFIGURATION THAT FITS YOUR WELLNESS GOALS.

The Vielight MIP supports brain and systemic wellness in three available configurations. Compare their applicators, wavelengths and benefits to find the option that best suits your routine.



| MIP 470–633–655–810                                                                      | MIP 633–810                                         | MIP 655–810                                         |
|------------------------------------------------------------------------------------------|-----------------------------------------------------|-----------------------------------------------------|
| 4 Applicators                                                                            | 2 Applicators                                       | 2 Applicators                                       |
| Wavelengths                                                                              | Wavelengths                                         | Wavelengths                                         |
| 810 nm (near-infrared)<br>633 nm<br>655 nm (laser)<br>470 nm                             | 810 nm (near-infrared)<br>633 nm                    | 810 nm (near-infrared)<br>655 nm (laser)            |
| Benefits                                                                                 | Benefits                                            | Benefits                                            |
| Gentle brain stimulation<br>Immunity<br>Oxygenation<br>Intranasal microbial sanitization | Gentle brain stimulation<br>Immunity<br>Oxygenation | Gentle brain stimulation<br>Immunity<br>Oxygenation |



# Designed for Everyday Use


## SIMPLE. PORTABLE. AUTOMATIC.

The Vielight MIP is designed to fit easily into an everyday wellness routine.



### 25-minute automatic sessions

The controller automatically ends each session after 25 minutes.



### Use one or two applicators

Connect the wavelength combination selected for your routine.



### Portable and battery operated

The lightweight controller operates with one AA alkaline or rechargeable NiMH battery.



### Automatic applicator detection

The controller recognizes the connected applicators and runs them accordingly.

